

# BLOOD MIXER (ROLLER MIXER)

“ISO 9001:2008 & ISO 13485:2003 ”  
CERTIFIED COMPANY  
“CE MARKED PRODUCTS”

This is an equipment to make blood homogeneous. It is shake-less and gentle mixing process. It is essential pre-process for perfect result of RBC, WBC and platelets counting on Hematology Analyzer.

## Model No : BM - 11

Roller Size : 4 Nos. X 230mm Lg.  
Mains Input : 230v, 50Hz / 110v, 60Hz, AC  
RPM : 30 RPM / 36 RPM  
Mixing Time : 15 to 20 Minutes  
Weight : 2.0 Kg. Net  
Capacity : Vacutainer 6 Nos.(of - 4ml)  
Or Bulb of 5ml x 15 Nos

