

# DIFFERENTIAL BLOOD CELL COUNTER

“ISO 9001:2008 & ISO 13485:2003”  
CERTIFIED COMPANY  
“CE MARKED PRODUCTS”

( With Computer Connectivity & Memory )

## Model No : DC - 33

Mains Input	: 230V/50Hz   110V/60Hz A.C. Supply
Body	: ABS Body
LCD Display	: 2-LINE 16 Character With Back Light
Weight	: 1 Kg.
Overall Size	: 155.0 (L) x 205 (W) x 91 (H)mm
Cell Counting Key	: 8 KEY
Total Key	: 12 KEY
Memory	: 50 Records
Computer	: Connectivity - USB port
Keypad	: Tectile switch with raised keytops



Computer Connectivity by “USB” port



LABLINE EQUIPMENTS PVT. LTD.  
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,  
BARODA 390010, GUJARAT, INDIA.  
PH. NO. +91 265 3927448 / 49 / 50  
FAX NO. +91 265 2632151  
E-MAIL : haresh@lablineindia.com

VISIT US AT :- [www.lablineindia.com](http://www.lablineindia.com)