

LABLINE'S DIGITAL CENTRIFUGES RANGE

MOTOR UNIVERSAL (AC/DC CARBON BRUSH TYPE)



RT-01



RT-02



RT-03



RT-04



RT-05



RT-06



RT-07



RT-08



RHC-12



RHC-24



MICRO TUBE



RD-11

DIGITAL CENTRIFUGE AND				HEMATOCRIT		CENTRIFUGE	
MODEL NO	TIMER	ROTOR	ROTOR NO	MAX RPM	SPEED REGULATOR	MAX RCF	NET WEIGHT
CF-600	YES	12 TUBE x 5 ML SWING OUT	RT-02	3800	STAPLES ELECTRONIC	2000	14.5 kg
	YES	8 TUBE x 15 ML ANGULAR	RT-03	3800	STAPLES ELECTRONIC	2000	14.5 kg
CF-608	YES	8 TUBE x 15 ML SWING OUT	RT-01	3200	STAPLES ELECTRONIC	1800	14.5 kg
CF-611	YES	8 TUBE x 15 ML ANGULAR	RT-03	5200	STAPLES ELECTRONIC	3300	14.5 kg
CF-622	YES	12 TUBE x 15 ML ANGULAR	RT-04	4900	STAPLES ELECTRONIC	3500	14.5 kg
CF-633	YES	16 TUBE x 15 ML ANGULAR	RT-05	4800	STAPLES ELECTRONIC	3500	14.5 kg
CF-644	YES	24 TUBE x 5 ML ANGULAR	RT-06	5500	STAPLES ELECTRONIC	2900	14.5 kg
CF-655	YES	16 TUBE x 15 ML SWING OUT	RT-07	3500	STAPLES ELECTRONIC	1900	14.5 kg
	YES	24 TUBE x 5 ML SWING OUT	RT-08	4000	STAPLES ELECTRONIC	1950	14.5 kg
HCF	YES	12 Capillaries of QBC - II Plus System	RHC-12	12000	STAPLES ELECTRONIC	15300	14.5 kg
	YES	24 Capillaries of size 75 mm x 1 mm bore	RHC-24	12000	STAPLES ELECTRONIC	15300	14.5 kg
	YES	24 TUBE x 1.5 ML 24 TUBE x 2.0 ML	MICRO TUBE	12000	STAPLES ELECTRONIC	13500	14.5 kg



LABLINE EQUIPMENTS PVT. LTD.
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,
BARODA 390010, GUJARAT, INDIA.
PH. NO. +91 265 3927448 / 49 / 50
FAX NO. +91 265 2632151
E-MAIL : haresh@lablineindia.com

VISIT US AT :- www.lablineindia.com