

# HEMATOCRIT CENTRIFUGE

(digital timer & RPM)

“ISO 9001:2008 & ISO 13485:2003”  
CERTIFIED COMPANY  
“CE MARKED PRODUCTS”



## Model No : HCF-12/24

|                 |  |
|-----------------|--|
| Capacity        | : 12 / 24 Nos. Capillary Tube          |
| Main input      | : 220v , 50Hz / 110v , 60Hz            |
| Max. Speed      | : More than 12000 R.P.M                |
| Timer           | : Digital Timer (0 to 60 minutes)      |
| Size            | : (W) 355 x (D) 410 x (H) 280 mm       |
| Weight          | : 15.0 Kg. Net                         |
| Motor           | : Universal (AC/DC Carbon Brush Type)  |
| RPM Indicator   | : Digital Up to 5 Digits : LCD Display |
| Speed Regulator | : Step less Electronic                 |

Note :- You can select any one rotor or all. The rotors are interchangeable in same device.



### RHC - 12

For 12 Capillaries of  
QBC - II Plus System  
R.P.M : More than 12000  
R.C.F : 15300g  
Material : Aluminum




### RHC - 24


For 24 Capillaries of size  
75 mm x 1 mm bore  
R.P.M : More than 12000  
R.C.F : 15300g  
Material : Aluminum



Micro Tube 1.5 ml / 2.0 ml  
R.P.M : More than 12000  
R.C.F : 13500g  
Material : Aluminum



Front View



Back View

## READING DEVICE RD - 11

Size : (L) 262 x (W) 110 x (H) 18  
Body : Aluminum & Plastic  
Weight : 0.25 Kg



LABLINE EQUIPMENTS PVT. LTD.  
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,  
BARODA 390013, GUJARAT, INDIA.  
PH. NO. +91 265 3927448 / 49 / 50  
FAX NO. +91 265 2632151  
E-MAIL : haresh@lablineindia.com

VISIT US AT :- [www.lablineindia.com](http://www.lablineindia.com)