

# CENTRIFUGE

“ISO 9001:2008 & ISO 13485:2003 ”  
 CERTIFIED COMPANY  
 “CE MARKED PRODUCTS”

## Model No : CF - 511    CF - 512

	4 x 15ml	6 x 15ml
Capacity	: 4 Nos. Tubes (Swing Out Head)	: 6 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz   110v,60Hz	: 220v,50Hz   110v,60Hz
Max. Speed	: More Than 3200 RPM	: More Than 3200 RPM
Max. R.C.F.	: 1600 x g	: 1700 x g
W x D x H	: 317 x 317 x 300mm	: 317 x 317x 300mm
Weight	: 6.0 Kg. Net	: 6.0 Kg. Net
Motor	: Universal ( AC/DC Carbon Brush type )	: Universal ( AC/DC Carbon Brush type )
Body	: Aluminum Body	: Aluminum Body



**Without Timer**



**Without Timer**

## Model No : CF - 513    CF - 514

	8 x 15ml	12 x 5ml
Capacity	: 8 Nos. Tubes (Swing Out Head)	: 12 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz   110v,60Hz	: 220v,50Hz   110v,60Hz
Max. Speed	: More Than 3200 RPM	: More Than 3200 RPM
Max. R.C.F.	: 1800 x g	: 1400 x g
W x D x H	: 317 x 317 x 300mm	: 317 x 317x 300mm
Weight	: 6.0 Kg. Net	: 6.0 Kg. Net
Motor	: Universal ( AC/DC Carbon Brush type )	: Universal ( AC/DC Carbon Brush type )
Body	: Aluminum Body	: Aluminum Body



LABLINE EQUIPMENTS PVT. LTD.  
 614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,  
 BARODA 390010, GUJARAT, INDIA.  
 PH. NO. +91 265 3927448 / 49 /50  
 FAX NO. +91 265 2632151  
 E-MAIL : haresh@lablineindia.com

VISIT US AT :- [www.lablineindia.com](http://www.lablineindia.com)

# CENTRIFUGE

“ISO 9001:2008 & ISO 13485:2003 ”  
CERTIFIED COMPANY  
“CE MARKED PRODUCTS”

## Model No : CF - 516    CF - 517

Capacity	<b>4 x 15ml</b> : 4 Nos. Tubes (Swing Out Head)	<b>6 x 15ml</b> : 6 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz   110v,60Hz	: 220v,50Hz   110v,60Hz
Max. Speed	: More Than 3200 RPM	: More Than 3200 RPM
Max. R.C.F.	: 1600 x g	: 1700 x g
W x D x H	: 317 x 317 x 300mm	: 317 x 317x 300mm
Weight	: 6.2 Kg. Net	: 6.2 Kg. Net
Motor	: Universal ( AC/DC Carbon Brush type )	: Universal ( AC/DC Carbon Brush type )
Body	: Aluminum Body	: Aluminum Body
Timer	: 0 to 30 Minutes	: 0 to 30 Minutes



**With Timer  
Auto Stop model**



**With Timer  
Auto Stop model**

## Model No : CF - 518    CF - 519

Capacity	<b>8 x 15ml</b> : 8 Nos. Tubes (Swing Out Head)	<b>12 x 5ml</b> : 12 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz   110v,60Hz	: 220v,50Hz   110v,60Hz
Max. Speed	: More Than 3200 RPM	: More Than 3200 RPM
Max. R.C.F.	: 1800 x g	: 1400 x g
W x D x H	: 317 x 317 x 300mm	: 317 x 317x 300mm
Weight	: 6.3 Kg. Net	: 6.3 Kg. Net
Motor	: Universal ( AC/DC Carbon Brush type )	: Universal ( AC/DC Carbon Brush type )
Body	: Aluminum Body	: Aluminum Body
Timer	: 0 to 30 Minutes	: 0 to 30 Minutes



LABLINE EQUIPMENTS PVT. LTD.  
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,  
BARODA 390010, GUJARAT, INDIA.  
PH. NO. +91 265 3927448 / 49 / 50  
FAX NO. +91 265 2632151  
E-MAIL : haresh@lablineindia.com

VISIT US AT :- [www.lablineindia.com](http://www.lablineindia.com)