

NEEDLE DESTROYER

“ISO 9001:2008 & ISO 13485:2003”
CERTIFIED COMPANY
“CE MARKED PRODUCTS”

This is an electrical equipment to destroy used NEEDLES by burning and destroy SYRINGE by cutting. We can protect nursing staff and society from transmission of diseases due to NEEDLE stick injuries.

Model No : ND - 22

Mains Input	: 230v, 50Hz / 110v, 60Hz , AC
Body	: Metal Body
Gauge Of Needles	: Can Burn 30 to 20 Gauge
Destruction Time	: 1.5 to 3 Seconds
Burning Temp.	: 1400° C to 1800° C
Weight	: 2.1 Kg. Net
Overall Size	: 155 (L) x 129 (W) x 86 (H) mm



WITH S. S. SYRINGE CUTTER



LABLINE EQUIPMENTS PVT. LTD.
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,
BARODA 390010, GUJARAT, INDIA.
PH. NO. +91 265 3927448 / 49 / 50
FAX NO. +91 265 2632151
E-MAIL : haresh@lablineindia.com

VISIT US AT :- www.lablineindia.com