

NEEDLE DESTROYER

“ISO 9001:2008 & ISO 13485:2003”
CERTIFIED COMPANY
“CE MARKED PRODUCTS”

This is an electrical equipment to destroy used NEEDLES by burning and destroy SYRINGE by cutting. We can protect nursing staff and society from transmission of diseases due to NEEDLE stick injuries.



ABS TOP COVER

Model No : ND - 55

| | |
|------------------|---------------------------------|
| Mains Input | : 230v, 50Hz 110v, 60Hz , AC |
| Body | : ABS Plastic + Metal |
| Gauge Of Needles | : Can Burn 30 to 20 Gauge |
| Destruction Time | : 1.5 to 3 Seconds |
| Burning Temp. | : 1400° C to 1800° C |
| Weight | : 1.9 Kg. Net |
| Overall Size | : 160 (L) x 133 (W) x 92 (H) mm |



LABLINE EQUIPMENTS PVT. LTD.
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,
BARODA 390010, GUJARAT, INDIA.
PH. NO. +91 265 3927448 / 49 / 50
FAX NO. +91 265 2632151
E-MAIL : haresh@lablineindia.com

VISIT US AT :- www.lablineindia.com