

# ROTARY MIXER

“ISO 9001:2008 & ISO 13485:2003”  
CERTIFIED COMPANY  
“CE MARKED PRODUCTS”



## Model No : RM - 11

Mains Input	: 230V/50Hz   110V/60Hz A.C. Supply
Body	: ABS + METAL Body
Motion	: Tilting - fix degree at 38°
Speed	: 15 RPM
Weight	: 1.2 Kg.
Size of Device	: 188 (L) x 221 (W) x 260 (H)mm
Operation	: Quiet - No noise
Standard configuration	: 4ml Vacutainer (EDTA) Tubes - 8 nos And ESR TUBE (8mm) - 8 nos. (Total-16 nos)
Special configuration	: 16 Tube of either type or mix as required.



LABLINE EQUIPMENTS PVT. LTD.  
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,  
BARODA 390010, GUJARAT, INDIA.  
PH. NO. +91 265 3927448 / 49 / 50  
FAX NO. +91 265 2632151  
E-MAIL : haresh@lablineindia.com

VISIT US AT :- [www.lablineindia.com](http://www.lablineindia.com)