

# ROCKING PLATFORM

“ISO 9001:2008 & ISO 13485:2003”  
CERTIFIED COMPANY  
“CE MARKED PRODUCTS”

Note:- We provide both types of rubber pads with device. Rubber pads are interchangeable.

Capacity : 24 Vaccutainer tube of 4 ml capacity  
Motion : Tilting - 15 Degree from horizontal  
on both side (Total 30 Degree)  
Main input : 220V / 50Hz Or 110V / 60 Hz  
Speed : 15 rocks / minute  
Platform size: 172 mm X 151 mm  
Weight : 2.2 Kg. Net  
Overall Size : 348 (L) x 198 (W) x 103 (H) mm

Model No : RP - 11



Model No : RP - 11



Capacity : 2-Tray of Western Blot  
Motion : Tilting - 15 Degree from horizontal  
On both side (Total 30 Degree)  
Main input : 220V / 50Hz Or 110V / 60 Hz  
Speed : 15 rocks / minute  
Platform size: 172 mm X 151 mm  
Weight : 2.2 Kg. Net  
Overall Size : 348 (L) x 198 (W) x 103 (H) mm



LABLINE EQUIPMENTS PVT. LTD.  
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,  
BARODA 390010, GUJARAT, INDIA.  
PH. NO. +91 265 3927448 / 49 / 50  
FAX NO. +91 265 2632151  
E-MAIL : haresh@lablineindia.com

VISIT US AT :- [www.lablineindia.com](http://www.lablineindia.com)