

CENTRIFUGE

“ISO 9001:2008 & ISO 13485:2003 ”
CERTIFIED COMPANY
“CE MARKED PRODUCTS”

Model No : CF - 11 CF - 22

Capacity	4 x 15ml : 4 Nos. Tubes (Swing Out Head)	6 x 15ml : 6 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz 110v,60Hz	: 220v,50Hz / 110v,60Hz
Max. Speed	: More Than 3100 RPM	: More Than 3100 RPM
Max. R.C.F.	: 1750 x g	: 1750 x g
W x D x H	: 300 x 300 x 245mm	: 300 x 300 x 245mm
Weight	: 4.0 Kg. Net	: 4.5 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)	: Universal (AC/DC Carbon Brush type)



Without Timer

Model No : CF - 33 CF - 44

Capacity	4 x 15ml : 4 Nos. Tubes (Swing Out Head)	6 x 15ml : 6 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz / 110v,60Hz	: 220v,50Hz / 110v,60Hz
Max. Speed	: More Than 3100 RPM	: More Than 3100 RPM
Max. R.C.F.	: 1750 x g	: 1750 x g
Timer	: 0 to 30 min.	: 0 to 30 min.
W x D x H	: 300 x 300 x 245mm	: 300 x 300 x 245mm
Weight	: 4.4 Kg. Net	: 4.5 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)	: Universal (AC/DC Carbon Brush type)



With Timer
Auto Stop model



LABLINE EQUIPMENTS PVT. LTD.
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,
BARODA 390010, GUJARAT, INDIA.
PH. NO. +91 265 3927448 / 49 / 50
FAX NO. +91 265 2632151
E-MAIL : haresh@lablineindia.com

VISIT US AT :- www.lablineindia.com