



“ISO 9001:2008 & ISO 13485:2003”  
 CERTIFIED COMPANY  
 “CE MARKED PRODUCTS”

# CENTRIFUGE

With Maintenance free - Brushless Motor  
 Microprocessor Controlled Digital Programmeble Timer



- Rotor : Swing out type [Horizontal]
- Capacity : 6 x 10 ml [17 Ø x 100 mm]
- Rubber cushion : For 6 x 5ml [EDTA / vaccute tube]
- Main input : 220v , 50Hz / 110v , 60Hz
- Speed : Fix - single speed machine  
 : 2850 RPM (220V/50Hz)  
 : 3400 RPM (110V/60Hz)
- Timer : Digital Timer (1 to 30 minutes)
- Size : (W) 325 x (D) 293 x (H) 230 mm
- Weight : 8 Kg. Net
- Motor : **Brush less motor**
- R.C.F : 1200g / 1600g
- Top cover : Transparent Poly Carbonat

## SAFETY FEATURES

1. Door open :- Device will not start
2. Door Latch :- Door remain locked till  
 Rotation Slow down or stop
3. Over current:- motor protected by 4 AMP MCB
4. Over heat :- Motor protected by Auto Thermal  
 Cutoff against over heat
5. Cycle counter:- Available with or without counter  
 It is digital electronics counter with LCD display

## Rotor with Holder



6 x 10 ml  
 swing out Type Rotor

MODEL	TYPE	MODEL	TYPE
Swing 60/CC	Programable Timer With Digital Cycle Counter	Swing 60 FT/CC	Fix Time Device With Digital Cycle Counter
Swing 60/WCC	Programable Timer Without Digital Cycle Counter	Swing 60 FT/WCC	Fix Time Device Without Digital Cycle Counter