

NEEDLE DESTROYER

This is an electrical equipment to destroy used NEEDLES by burning and destroy SYRINGE by cutting. We can protect nursing staff and society from transmission of diseases due to NEEDLE stick injuries.

Model No : ND - 22

Mains Input	: 230v, 50Hz or 110v, 60Hz , AC
Body	: Metal Body
Gauge Of Needles	: Can Burn 30 to 20 Gauge
Destruction Time	: 1.5 to 3 Seconds
Burning Temp.	: 1400° C to 1800° C
Weight	: 2.1 Kg. Net
Overall Size	: 155 (L) x 129 (W) x 86 (H) mm
Cutter	: S.S. Blade



Metalic Top Cover



ABS TOP COVER

Model No : ND - 55

Mains Input	: 230v, 50Hz or 110v, 60Hz , AC
Body	: ABS Plastic + Metal
Gauge Of Needles	: Can Burn 30 to 20 Gauge
Destruction Time	: 1.5 to 3 Seconds
Burning Temp.	: 1400° C to 1800° C
Weight	: 1.9 Kg. Net
Overall Size	: 160 (L) x 133 (W) x 92 (H) mm
Cutter	: S.S. Blade

Model No : ND - 66

Mains Input	: 230v, 50Hz or 110v, 60Hz
Body	: Metal Body
Gauge Of Needles	: Can Burn 30 to 20 Gauge
Destruction Time	: 1.5 to 3 Seconds
Burning Temp.	: 1400° C to 1800° C
Weight	: 1.9 Kg. Net
Overall Size	: 155 x 128 x 80 mm
Cutter	: S.S. Blade

