

## DOCTOR CENTRIFUGE

### MOTOR UNIVERSAL (AC/DC CARBON BRUSH TYPE)

**A**



Model	: CF - D11
Capacity	: 8 x 15 ml Fixed Angled Rotor
Main input	: 220v , 50Hz / 110v , 60Hz
Speed	: 4000 R.P.M-Maximum
Speed Regulator	: Electronics Steples Speed Regulator
Size	: (W) 300 x (D) 300 x (H) 320 mm [PACKING]
Weight	: 6.5 Kg. Net
Motor	: Universal (AC/DC carbon brush type)
R.C.F	: 1550g

**B**

Model	: CF - D22
Capacity	: 8 x 15 ml Fixed Angled Rotor
Main input	: 220v , 50Hz / 110v , 60Hz
Speed	: 4000 R.P.M-Maximum
Timer	: CF-D11 - With Out Timer : CF-D22 - Digital Timer (1 to 30 minutes)
Speed Regulator	: Electronics Steples Speed Regulator
Size	: (W) 300 x (D) 300 x (H) 320 mm [PACKING]
Weight	: 6.5 Kg. Net
Motor	: Universal (AC/DC carbon brush type)
R.C.F	: 1550g

