

BLOOD MIXER (ROLLER MIXER)

This is an equipment to make blood homogeneous. It is shake-less and gentle mixing process. It is essential pre-process for perfect result of RBC, WBC and platelets counting on Hematology Analyzer.

Model No : BM - 11

- Roller Size : 4 Nos. X 230mm Lg.
- Mains Input : 230v, 50Hz or 110v, 60Hz, AC
- RPM : 30 RPM
- Mixing Time : 15 to 20 Minutes
- Weight : 2.0 Kg. Net
- Capacity : Vacutainer 6 Nos.(of - 4ml)
Or Bulb of 5ml x 15 Nos



Model No : BM - 22

- Roller Size : 4 Nos. X 340mm Lg.
- Mains Input : 230v, 50Hz or 110v, 60Hz, AC
- RPM : 30 RPM
- Mixing Time : 15 to 20 Minutes
- Weight : 2.4 Kg. Net
- Capacity : Vacutainer 12 Nos.(of - 4ml)
Or bulb of 5ml x 24 Nos