

MOTOR UNIVERSAL (AC/DC CARBON BRUSH TYPE)

A



B



A

SMALL CENTRIFUGE MACHINE

MODEL NO	TIMER	ROTOR	MAX RPM	SPEED REGULATOR	MAX RCF	WEIGHT
CF-11	NO	4 TUBE x 15 ML SWING OUT	4000	STEPLESS ELECTRONICS	2275	4.0 kg
CF-22	NO	6 TUBE x 15 ML SWING OUT	3400	STEPLESS ELECTRONICS	1800	4.5 kg
CF-33	YES	4 TUBE x 15 ML SWING OUT	4000	STEPLESS ELECTRONICS	2275	4.3 kg
CF-44	YES	6 TUBE x 15 ML SWING OUT	3400	STEPLESS ELECTRONICS	1800	4.8 kg

B

MEDIUM CENTRIFUGE MACHINE

MODEL NO	TIMER	ROTOR	ROTOR NO	MAX RPM	SPEED REGULATOR	MAX RCF	WEIGHT
CF-511	NO	4 TUBE x 15 ML SWING OUT	—	4000	STEPLESS ELECTRONICS	2275	6.0 kg
CF-512	NO	6 TUBE x 15 ML SWING OUT	—	3400	STEPLESS ELECTRONICS	1800	6.0 kg
CF-513	NO	8 TUBE x 15 ML SWING OUT	RT-01	3400	STEPLESS ELECTRONICS	1875	6.0 kg
CF-514	NO	12 TUBE x 5 ML SWING OUT	RT-02	4000	STEPLESS ELECTRONICS	2000	6.0 kg
CF-516	YES	4 TUBE x 15 ML SWING OUT	—	4000	STEPLESS ELECTRONICS	2275	6.3 kg
CF-517	YES	6 TUBE x 15 ML SWING OUT	—	3400	STEPLESS ELECTRONICS	1800	6.3 kg
CF-518	YES	8 TUBE x 15 ML SWING OUT	RT-01	3400	STEPLESS ELECTRONICS	1875	6.3 kg
CF-519	YES	12 TUBE x 5 ML SWING OUT	RT-02	4000	STEPLESS ELECTRONICS	2000	6.3 kg