

ROTOR FOR DIGITAL CENTRIFUGES CF - 600 / SC - 888 / CF - 800 SERIES

RT-01



RT-02



RT-03



RT-04



RT-05



RT-06



RT-07



RT-08



RT-09



RT-10



RT-11



RT-12



RT-13



RHC-12




RHC-24



MICRO TUBE



RD-11



SR. NO	ROTOR NO.	MODEL NO			TYPE OF ROTOR	MAX RPM	MAX RCF
		CF - 600 SERIES	SC - 888 SERIES	CF - 800 SERIES			
01	RT - 01	CF - 608	SC - 888	CF - 811	8 TUBE x 15 ML SWING OUT	3500	1850
02	RT - 02	CF - 600	SC - 888	CF - 822	12 TUBE x 05 ML SWING OUT	4000	2000
03	RT - 03	CF - 611	SC - 888	CF - 833	8 TUBE x 15 ML ANGULAR	5000	3150
04	RT - 04	CF - 622	SC - 888	CF - 833	12 TUBE x 15 ML ANGULAR	5000	3450
05	RT - 05	CF - 622	SC - 888	CF - 833	16 TUBE x 15 ML ANGULAR	5000	3724
06	RT - 06	CF - 611	SC - 888	CF - 833	24 TUBE x 05 ML ANGULAR	5000	2900
07	RT - 07	CF - 655	SC - 888	CF - 811	16 TUBE x 15 ML SWING OUT	3500	1975
08	RT - 08	CF - 655	SC - 888	CF - 822	24 TUBE x 05 ML SWING OUT	4000	2024
09	RT - 09	CF - 600	SC - 888	CF - 811	12 TUBE x 07 ML SWING OUT	3500	1975
10	RT - 10	CF - 600	SC - 888	CF - 811	2 TUBE x 50 ML SWING OUT	3500	1975
11	RT - 11	CF - 655	SC - 888	CF - 811	24 TUBE x 07 ML SWING OUT	3500	1975
12	RT - 12	CF - 655	SC - 888	CF - 811	4 TUBE x 50 ML SWING OUT	3500	1975
13	RT - 13	NA	SC - 888	CF - 822	24 TUBE x 15 ML ANGULAR	4000	2437
14	RT - 14	NA	SC - 888	CF - 822	16 TUBE x 15 ML SWING OUT (Metallic)	4500	3150
15	RHC - 12	HCF	NA	CF - 844	12 CAPILLARIES OF QBC - II PLUS SYSTEM	12000	15300
16	RHC - 24	HCF	NA	CF - 844	24 CAPILLARIES OF SIZE 75 MM x 1 MM BORE	12000	15300
17	MICRO TUBE	HCF	NA	CF - 844	24 TUBE x 1.5 ML / 2.0 ML ANGULAR	12000	15300