

DRY BATH INCUBATOR

Model No : DBI - 11 / DBI - 66



Mains Input	: 230V / 50 Hz OR 110V / 60Hz A.C. Supply
W x D x H	: 168 (W) x 232 (D) x 93 (H)
Body	: ABS Plastic Body & Metal Base (DBI-11) : Metal Body (DBI-66)
Weight	: 2.0 Kg. Net
Temp.Controller	: Digital Type
Temp. Range	: Ambient +5°C To 70°C (DBI-11) : Ambient +5°C To 110°C (DBI-66)
Accuracy	: $\pm 0.5^{\circ}\text{C}$ @ 37°C, $\pm 3\%$ Except 37°C
Timer	: 1 To 999 (Digital Type)
Option Temp.	: 85°C for DBI-11
Display	: LCD TYPE

Block - 1




12.7 MM 24 HOLE

Block - 2




10.5 MM 28 HOLE

Block - 3




**18.5 MM - 4 HOLE
12.7 MM - 6 HOLE
15.5 MM - 5 HOLE
10.5 MM - 7 HOLE**

Block - 4




**12.7 X 12.7 MM
20 SQUAR HOLE**

Block - 5




**1.5 ML - 24 HOLE
TAPER HOLE 11.0 DIA**

Block - 6




**16.5MM DIA
10 HOLE**

Block - 7



**17.5MM DIA
10 HOLE**

Block - 8



**0.1 ML / 0.2 ML
36 HOLE**