

With Maintenance free - Brushless Motor
Microprocessor Controlled Digital Programmable **Timer & Variable Speed.**



MODEL :- SC -101

Rotor : As Per Selected Rotor
Main input : 220v, 50Hz or 110v, 60Hz
Speed : **Variable Up To 4000 RPM**
Timer : Digital Timer (1 to 30 minutes)
Size : (W) 325 x (D) 293 x (H) 230 mm
Weight : 8 Kg. Net
Motor : **BRUSH LESS MOTOR**
Top cover : Almost Unbreakable polycarbonet
Digital display : LCD display of RPM , RCF , TIME

SAFETY FEATURES

1. The device will not start if door is open
2. The device will shut down if door is open during processing. Do not open the door if rotor in motion
3. Motor protected by Auto Thermal Cutoff against overheat

FOR ABOVE MACHINE PRICE INCLUDE SWING OUT ROTOR + 2 NOS. ADAPTOR FOR 8 x 15 ML + 2 NOS ADAPTOR FOR 12 x 05 ML. **OR 8 x 15 ML ANGULAR HEAD ROTOR, **ANY ONE.** FOR SWING OUT TYPE ROTOR 12 x 7 ML & 2 x 50 ML ADAPTOR CAN PURCHASE EXTRA.**

ROTORS FOR BOTH MACHINE

RT-01



8 x 15ML SWING OUT
R.C.F : 2350g

RT-02



12 x 05ML SWING OUT
R.C.F : 1950g

RT-03




08 x 15ML ANGULAR HEAD
R.C.F : 2050g

RT-09



12 x 07ML SWING OUT
R.C.F : 2300g

RT-10



02 x 50ML SWING OUT
R.C.F : 2350g