

NON DIGITAL BRUSH LESS CENTRIFUGE

With Maintenance free - Brushless Motor

Model : CF - 521



WITH OUT TIMER

Rotor : Swing out type [Horizontal]
Capacity : As Per Selected Rotor
Main input : 220v, 50Hz or 110v, 60Hz
Speed : Variable Up To 4000 RPM
Size : (W) 320 x (D) 320 x (H) 300 mm
Weight : 8 Kg. Net
Motor : BRUSH LESS MOTOR

Model : CF - 522



WITH TIMER AUTO STOP

Rotor : Swing out type [Horizontal]
Capacity : As Per Selected Rotor
Main input : 220v, 50Hz or 110v, 60Hz
Speed : Variable Up To 4000 RPM
Timer : Digital Timer (1 to 30 minutes)
Size : (W) 320 x (D) 320 x (H) 300 mm
Weight : 8 Kg. Net
Motor : BRUSH LESS MOTOR

FOR ABOVE MACHINE PRICE INCLUDE SWING OUT ROTOR + 2 NOS. ADAPTOR FOR 8 x 15 ML + 2 NOS ADAPTOR FOR 12 x 05 ML. **OR** 8 x 15 ML ANGULAR HEAD ROTOR, **ANY ONE**.
FOR SWING OUT TYPE ROTOR 12 x 7 ML & 2 x 50 ML ADAPTOR CAN PURCHASE EXTRA.

ROTORS FOR BOTH MACHINE

RT-01



8 x 15ML SWING OUT
R.C.F : 2350g

RT-02



12 x 05ML SWING OUT
R.C.F : 1950g

RT-03




08 x 15ML ANGULAR HEAD
R.C.F : 2050g

RT-09



12 x 07ML SWING OUT
R.C.F : 2300g

RT-10



02 x 50ML SWING OUT
R.C.F : 2350g