



“ISO 9001:2015 & EN ISO 13485:2016”
CERTIFIED COMPANY
“CE MARKED PRODUCTS”

DRY BATH INCUBATOR

[FIX TEMPERATURE TYPE]

Model No : DBI - 22



Mains Input	: 230V/50 Hz OR 110V/60Hz A.C. Supply
W x D x H	: 195 (W) x 86 (D) x 80 (H)
Body	: ABS & Metal Body
Weight	: 1.5 Kg. Net
Temp. Controller	: Digitally Controlled
Temp. Fix	: 37 °C Std or 37 °C to 50 °C on request
Accuracy	: ± 0.5°C
Heating Block - 1	: 16 Holes Of 12.7 MM Dia
Heating Block - 2	: 16 Holes Of 1.5 ml Tube
Heating Block - 3	: 16 Holes Of 10.5 mm Dia

Model No : DBI - 44

Mains Input	: 230V/50Hz 110V/60Hz A.C. Supply
W x D x H	: 114 (w) x 90 (D) x 80 (H) MM
Body	: ABS & Metal Body
Weight	: 1.0 KG. NET
Temp. Controller	: Digitally Controlled
Temp. Fix	: 37 °C - Standard or 37 C to 50 °C on request
Accuracy	: ± 0.5°C
Heating Block - 1	: 10 Holes Of 1.5 ml Tube
Heating Block - 2	: 10 Holes Of 10.5 mm dia

