

With Maintenance free - Brushless Motor
Microprocessor Controlled Digital Programmable **Timer & Variable Speed.**



MODEL :- MC - 11 / MC - 22

Rotor : 8 X 1.5 ml / 8 X 2 ml /
8-2 line X 0.2 ml PCR strip

Main input : 220v, 50Hz or 110v, 60Hz

Speed : **Variable Up To 6000 RPM**

Timer : Digital Timer (1 to 60 minutes)

Size : (W) 172 x (D) 190 x (H) 145

Weight : 3 Kg Net.

Motor : BRUSH TYPE DC MOTOR

Top cover : Almost Unbreakable polycarbonet

Digital display : LCD display of RPM , TIME

SAFETY FEATURES

1. The device will not start if door is open
2. The device will shut down if door is open during processing. Do not open the door if rotor in motion