

# VORTEX MIXER

**Model No : VM - 11**

Mains Input : 230V/50Hz or 110V/60Hz  
 A.C. Supply

W x D x H : 160 (W) x 190 (D) x 115 (H)

Body : ABS Plastic Body & Metal Base

Weight : 3.0 Kg. Net

Speed : 300- 2900 RPM (Variable)



Operation modes : Touch Or Continuous

# DIFFERENTIAL BLOOD CELL COUNTER



**Model No : DC - 22**

Mains Input : 230V/50Hz or 110V/60Hz A.C. Supply

Body : ABS Body ( Light Weight )

LCD Display : 16 Character With Back Light

Weight : 1.0 Kg.

Overall Size : 155.0 (L) x 205 (W) x 91 (H)mm

Cell Counting Key : 8 KEY

Total Key : 12 KEY

Keypad : Tactile switch with raised keytops

Timer & Stop Upto 999 Minutes

( With Computer Connectivity & Memory )

**Model No : DC - 33**

Mains Input : 230V/50Hz or 110V/60Hz A.C. Supply

Body : ABS Body

LCD Display : 2-LINE 16 Character With Back Light

Weight : 1 Kg.

Overall Size : 155.0 (L) x 205 (W) x 91 (H)mm

Cell Counting Key : 8 KEY

Total Key : 12 KEY

Memory : 50 Records

Computer : Connectivity - USB port

Keypad : Tactile switch with raised keytops



Computer Connectivity by “USB” port