

DRY BATH INCUBATOR

[FIX TEMPERATURE TYPE]

Model No : DBI - 55 DUAL BLOCK



Mains Input	: 230V/50Hz 110V/60Hz A.C. Supply
W x D x H	: 130 (w) x 160 (D) x 92 (H) MM
Body	: ABS & Metal Body
Weight	: 1.3 kG. NET
Temp. Controller	: Digitally Controlled
Accuracy	: $\pm 0.5^{\circ}\text{C}$
Temp. Fix	: 30°C - FOR CASSET (ROOM TEMP REQUIRED BELOW 25°C) AND 50°C FOR MICRO TUBE
Heating Block - 1	: 10 Holes Of 1.5 ml Tube OR 10 Holes Of 2.0 ml Tube
Heating Block - 2	: 9 SLOTS OF 15(w) x 9(h) X 63(d) MM